

ORIPLAST®



Surgical suction

RELIABLE CLARIFICATION


Surgical suction

SURGICAL FLEXIBLE SUCTION SYSTEMS

Systems for every need


Surgical suction systems are designed for the suction of liquids and tissue debris. They are mainly used during operations to keep the surgical field clean and to maintain an overview or to remove the rinsing fluid that has been used.

Surgical suction systems remove solid and liquid substances quickly and reliably without clogging.



Our connection tubing s offer numerous connection options with a **large variety of Approaches**. The following approaches can be found in our range: Sherwood, funnel, silicone and fingertip.

Approach



Suction tubing

The length of the suction tubing s varies from 1 to 3 m and enables an optimal and reliable connection option with a large number of Approach s. They are also available in different lumen sizes. The transparent disposable tubing s are made of **high-quality, kink-resistant material**.

Fingertip

The fingertip can be used as a connection option. This is used in a regulating manner and, with the integrated stepped connector, offers further connection options.

Suction tubing

Surgical suction tubing s are suitable for the suction of body fluids. These are available in various designs and sizes.



Suction tubing

- Different approaches for maximum compatibility
- High-quality materials

eClass 34-24-13-05



STERILE EO

Item No.	PZN	Lenght (mm)	FG	Approach 1	Approach 2	Protective cap	DEHP	PVC	PU
VK00005	17483273	1800	18	Silicone 5 x 1,5 mm	Silicone 7 x 1,5 mm	-			100
VK00006	17483296	1800	18	Silicone	LL pos	1			50
VK00112	17484166	2100	25	Funnel	Funnel	-			50
VK00113	17484172	3000	25	Funnel	Funnel	-			50
VK00128	17484249	1600	25	Silicone	Silicone	-			50
VK00131	17484278	1600	25	Silicone	Fingertip	-			50
VK00516	17486308	2000	25	Sherwood	Sherwood	-			50
VK00519	17486320	3000	25	Sherwood	Sherwood	-			50
VK00523	17486337	3000	25	Sherwood	Sherwood	-			50
VK00524	17486343	3000	25	Sherwood	Silicone	-			50
VK00527	17486366	1600	25	Sherwood	-	-			50
VK00529	17486372	1600	25	Sherwood	Fingertip	-			50
VK00007	17483304	1000	27	Silicone	Silicone	-			50
VK00127	17484232	1300	27	Silicone	Silicone	-			50
VK00001	17483267	1800	33	Silicone	Fingertip (large)	2			50
VK00111	17484143	1900	40	Step connector	-	-			25



Connecting tubing

- Different approaches for maximum compatibility
- High-quality materials

eClass 34-24-13-05



STERILE EO

Item No.	PZN	Lenght (mm)	FG	Approach	DEHP	PVC	PU
VK00115	17484189	3000	25	-			50
VK01060	17489413	2200	25	-			50
VK00916	17488537	280	28	Y connector			50

Fingertip

The fingertip is an essential tool when suctioning secretions. The fingertip can be opened or closed with the fingertip, which enables suction regulation and thus gentle suction.



Fingertip

- with a round support
- with sealing plug

eClass 34-24-13-07



STERILE | EO

Item No.	PZN	Approach 1	Approach 2	DEHP	PVC	PU
VK00950	17488796	Step connector	Step connector			100
VK00952	17488810	Step connector	Luer plug			100

Aspiration

The Yankauer suction device is a tool for a gentle and efficient suction method without aspiration of the surrounding tissue. This can be connected to a fingertip, which allows the vacuum to be regulated by thumb or finger pressure.



Yankauer

- Length: 280mm
- Angled suction tip
- Without suction tubing

eClass 34-31-05-01



STERILE | EO

Item No.	PZN	FG	Size	DEHP	PVC	PU
VK00226	17484953	22	Large			100
VK00227	17484976	18	Medium			100
VK00228	17484982	12	Small			100

Adapter



Tubing adapter

- for connecting and adapting
- Sherwood approach on both sides

eClass 34-24-13-92



STERILE | EO

Item No.	PZN	Lenght (mm)	FG	DEHP	PVC	PU
VK00515	17486283	150	25			100