

ORIPLAST®



Pleural puncture

HIT PRECISE TARGETS

Pleural puncture

RECONSTRUCTION OF THE RESPIRATORY SYSTEM

Special puncture method

Utensils for pleural puncture are essential, especially in emergency medicine: They are used to immediately drain fluid from the chest (pleural effusion), for example, which impedes the patient's breathing. In addition to such a therapeutic thoracocentesis, the diagnostic pleural puncture is important, which serves to collect and analyze the effusion in a sterile manner.

Possible causes of a pleural effusion are, for example, inflammatory changes, tumor diseases, traumatic injuries and heart diseases. The puncture is performed under local anesthetic with the patient sitting or lying down.

The clear three-piece syringe is designed for quick and easy connection to the cannula with a **luer lock fitting**. The syringe has a volume of 50 ml.

Syringe


A 14G cannula is included in our pleural puncture set. This can be connected quickly and effortlessly thanks to the Luer connection and is immediately ready for use.

Cannula

The hub is suitable for a **secure connection** of a cannula.

Luer lock connector





Connection tubing

The transparent connection tubing has a Luer lock and offers numerous connection options for the secretion bag, which is used for the **hygienic and safe drainage** of bodily fluids.

Three-way tap extension

Easy to use thanks to a three-way tap that cannot be mixed up.

Secretion bag

You will find our secretion bags in volumes of 600 ml, 1500 ml and 2000 ml. All secretion bags are provided with an easy-to-read scale so that you can keep an eye on the volume of the drained liquid at all times.

Pleural puncture set

The compact set is used in the diagnostic and therapeutic puncture of the pleural space between the ribs and the lungs, e.g. in the presence of a pleural effusion.



Pleural puncture set

- Tubing length 900 mm
- NRV
- Blue three-way cock
- 50 ml Luer lock syringe
- Puncture cannula 14G
- Protective cap

eClass 34-29-90-90



STERILE | EO

Item No.	PZN	Volume (ml)	Extension (mm)	Bottom drain	Union nut	DEHP	PVC	PU
VK00134	17484284	1500	400	-	-			50
VK00136	17484309	2000	400	•	-			50
VK00140	17484315	1500	300	-	•			50

Accessories



Pleural puncture set

- Tubing length 900 mm
- Blue three-way cock
- 50 ml Luer lock syringe
- Protective cap

eClass 34-29-02-90



STERILE | EO

Item No.	PZN	Volume (ml)	NRV	Extension (mm)	Approach	Packaging	DEHP	PVC	PU
VK00135	17484290	1500	•	400	LL pos				50
VK01061	17489436	-	-	300	LL neg				20



Secretion exchange bag

- Tubing length 900 mm
- LL neg

eClass 34-29-03-92



STERILE | EO

Item No.	PZN	Volume (ml)	NRV	DEHP	PVC	PU
VK00918	17488543	1500	•			100
VK01070	17489502	1500	-			100